

# Utilisation de l'électrophorèse capillaire en mode frontal continu (FACCE) pour la caractérisation des interactions entre facteurs de transcription et séquences d'ADN

## cibles, un pas vers le développement d'un microsystème analytique.

The objective of this project is to address the detailed study of interactions between FTs and their DNA targets by developing a method to characterise these interactions. In this context, the FACCE approach (continuous front-end capillary electrophoresis) offers a very promising approach. Indeed, FACCE is a free medium electrophoretic technique (absence of gel) allowing the quantification of free ligand in pre-equilibrated substrate-ligand mixtures, and has the advantage of using small sample volumes (a few nL). The quantification of the free ligand, in mixtures with varying ligand concentrations, allows the establishment of an interaction isotherm that allows the determination of the interaction parameters between the ligand and the substrate, notably the stoichiometry and the association constant. This technique has already been successfully applied to the study of protein/protein or biopolymer/protein interactions for which affinity constants have been measured

## RESULTATS

We studied the interaction between the transcription factor ILR3/bHLH105, recently identified (by the BPMP partner) as playing a central role in the control of iron homeostasis in *Arabidopsis thaliana* (a plant model organism), and one of its target sequences. This interaction was previously validated in vivo (single hybrid in yeast) and in vitro by gel delay (EMSA: electromobility shift assay). Analyses conducted by capillary electrophoresis in continuous frontal mode, at physiological pH, made it possible to determine the interaction parameters between the ILR3 protein and the target DNA. A dissociation constant  $k_d = 2.49E-7$  M<sup>-1</sup> and a complex stoichiometry  $n = 1.3$  were found. Although there is a large uncertainty in the value of  $n$  due to the difficulty of determining the exact concentration of pure protein, the value obtained suggests a 1:1 complex in agreement with what is expected.

**Responsable :**

**Date de démarrage :** 01/01/2017

**Date de clôture :** 30/06/2018

**Montant :**

